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MOLECULAR SIGNATURE TEST

prism/RA

HEALTHCARE PROVIDER	CLINIC	PATIENT
Name: Robert Jones Requisition No.: 00000022 Collection Date: 00-00-2022 Collection Time: 00:00 PM Specimen Received: 00-00-2022 Report Date: 00-00-2022	Clinic Name: Anytown Community Rheumatology Clinic Address: 1234 Mainstreet, Suite A1, Anytown, PA, 00000 Phone: 000-000-0000 Fax: 000-000-0000	Name: Jane Doe DOB: 00/00/1950 Sex: Female Weight: 150 Height: 5' 7" Patient Global Assessment: 65

**PrismRA RESULT**

<p><b>PrismRA TNFi Result:</b></p> <p style="text-align: center;"><b>TNFi Inadequate Response Signature Not Detected</b></p> <p style="text-align: center;">This patient has a PrismRA Score of <b>9.2</b></p>	<p><b>Interpretive Criteria:</b></p> <div style="text-align: center; border: 1px solid #2e8b57; padding: 5px; width: fit-content; margin: 0 auto;"> <span style="font-size: 2em; color: #2e8b57;">✓</span> </div> <p style="text-align: center;"><b>Patient may be likely to respond to a TNFi</b></p> <p style="text-align: center;">&lt;<b>10.6</b>: Patient does not have a detectable molecular signature of inadequate response to TNFi therapies</p>
--	--

**PrismRA SCORE INTERPRETATION**

A patient with a PrismRA Score <10.6 does not have a detectable molecular signature of inadequate response to TNFi therapies. The PrismRA result is reported on a continuous 1-25 scale. The higher the score, the more likely the patient will have an inadequate response to TNFi therapies; the lower the score, the less likely the patient will have an inadequate response to TNFi therapies.

**TEST DESCRIPTION**

PrismRA is a molecular signature response classifier (MSRC) that uses gene expression features, clinical features and anti-cyclic citrullinated protein (anti-CCP) antibody to detect a signature of non-response to tumor necrosis factor- $\alpha$  inhibitors (TNFi) for patients with rheumatoid arthritis. PrismRA predicts the patient's likelihood of inadequately responding to all TNFi therapies. Total RNA was extracted from the blood sample and sequenced at Ambry Genetics using validated methods. Anti-CCP antibody measurements were performed at Scipher using validated methods. Final algorithmic analysis was performed at Scipher with PrismRA pipeline v1.1. Response is defined as achieving ACR50 at 6 months. The PrismRA result is reported on a continuous 1 to 25 scale. The higher the score, the more likely the patient will have an inadequate response to TNFi therapies; the lower the score, the less likely the patient will have an inadequate response to TNFi therapies. However, a low score does not ensure a positive response to TNFi therapies.

**COMMENTS**

The PrismRA result is intended for informational purposes only and does not constitute a recommendation. Medical decisions should be made by the patient's provider. This service is developed and its performance characteristics determined by Scipher Medicine. This test has not been certified or approved by the US Food and Drug Administration. The laboratory is certified under the Clinical Laboratory Improvement Amendments (CLIA) as qualified to perform high complexity clinical testing for all test specifications and related publications please visit www.prismra.com

**Scipher Medicine**  
8000 Glenview, Suite 200  
Darien, NC 27008  
855-241-3127 www.prismra.com

**Ambry Genetics**  
Apo Valley, CA 92655  
Laboratory CLIA No. 000808484

**Laboratory CLIA Number:**  
240309014  
**Laboratory Director:**  
Dawn Gattuso, PhD  
Scott Rappaport, PhD, DABCCP

©2022 Scipher Medicine Corporation. All rights reserved. Scipher Medicine and PrismRA are trademarks of Scipher Medicine Corporation. All other trademarks are the property of their respective owners.

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
1. Arthritis Foundation: Rheumatoid Arthritis and Heart Disease. [www.arthritis.org/health-wellness/about-arthritis/related-conditions/other-diseases/rheumatoid-arthritis-heart-disease](http://www.arthritis.org/health-wellness/about-arthritis/related-conditions/other-diseases/rheumatoid-arthritis-heart-disease). Accessed 4/22/22. 2. Jin, Y., Desai, R.J., Liu, J. et al. Factors associated with initial or subsequent choice of biologic disease-modifying antirheumatic drugs for treatment of rheumatoid arthritis. *Arthritis Res Ther* 19, 159 (2017). 3. Mellors T, et al. Clinical validation of a blood-based predictive test for stratification of response to tumor necrosis factor inhibitor therapies in rheumatoid arthritis patients. *Network and Systems Medicine* 2020; 3(1):1–14. 4. Cohen, S., Wells, A.F., Curtis, J.R. et al. A Molecular Signature Response Classifier to Predict Inadequate Response to Tumor Necrosis Factor- $\alpha$  Inhibitors: The NETWORK-004 Prospective Observational Study. *Rheumatol Ther* 8, 1159–1176 (2021). 5. Creaky Joints: RA Testing Patient Guideline. [https://creakyjoints.org/www/wp-content/uploads/CreakyJoints\\_RA\\_Testing\\_Patient\\_Guideline.pdf](https://creakyjoints.org/www/wp-content/uploads/CreakyJoints_RA_Testing_Patient_Guideline.pdf). Accessed 6/1/22. 6. Strand V, Cohen SB, Curtis JR, et al. Clinical utility of therapy selection informed by predicted nonresponse to tumor necrosis factor- $\alpha$  inhibitors: an analysis from the Study to Accelerate Information of Molecular Signatures (AIMS) in rheumatoid arthritis. *Expert Rev Mol Diagn.* 2022;22(1):101-109. 7. Fraenkel, L., et al. (2021), American College of Rheumatology Guideline for the Treatment of Rheumatoid Arthritis. *Arthritis Care Res*, 73: 924-939.